



Dr. _____ Delivery Date _____

Mr./Ms. _____


Abutment Material Choice

- Titanium Abutment
- Zirconia w/Ti-Insert
- All-Zirconia Abutment
- _____

Parallel Abutment

- NO YES (indicate which abutments will have restorations splinted together for insertion)

Abutment Margin Depth



*If left blank, default values will be used

- Buccal/facial:** 1.0 mm subgingival
- Distal:** 0.75 mm subgingival
- Mesial:** 0.75 mm subgingival
- Lingual:** 0.5 mm subgingival

Abutment Margin Design

- Shoulder for all-ceramic
- Chamfer for PFM/BruxZir

Abutment Emergence Profile

- Surgical Placement
- Tissue Displacement
- No Tissue Displacement

Additional Instructions:

Select Abutment Type

NobelReplace

- NP WP
- RP Ø6

Neoss

- 4.0 mm

Branemark

- RP

NobelActive

- NP RP

Astra

- 3.0 4.5/5.0
- 3.5/4.0

Biomet 3i

- 3.4 mm 5.0 mm
- 4.1 mm 6.0 mm

PrimaConnex

- SD WD
- RD

Straumann Bone Level

- NC RC

Zimmer Screw-Vent

- 3.5 mm 5.7 mm
- 4.5 mm